

SightRPI

Fall Detection System

User Manual

You are 6 steps away from a safe life!

How System Works – At home

The solution that cares all about your safety and convenience.




Let us install for you!

Please cooperate with our work staff to install and test system embedded inside wall and ceiling.


How System Works – On Your Phone

The solution that cares all about your safety and convenience.




Download our App from app store.

- If you use **Android** phone, you can download it from **Google Play** store or from our website.
- If you use **iOS** devices, you can download it from **App Store**.



Connect Bluetooth.

- Turn on Bluetooth on your phone and click on “Connect” button circled by yellow rectangle.
- App will automatically search for Bluetooth available for connection. Choose our system Bluetooth “SightRP1”.
- “Connect” will change to red “Disconnect” (as circled in red circle). If not, please check Bluetooth setting or contact us.



SightRP1


4

If you are user's
emergent contact,
skip to **Step 6**.

On user's phone,

Add emergent contact.


- Type in your emergent contact phone number in the slot circled by yellow rectangle.
- Press "SET".



5

Start safe life!

- Click on user mode "USER(LIGHT CONTROL)".
- To adjust light color
 - Adjust Red, Green, Blue attributes.
 - Click "change" right below.
- To switch ON/OFF light
 - Swipe "LIGHT" button.
 - Click "SEND" below.



SightRPI

6

On emergent contact's phone,

Free Mode

- Our system works correctly and there is no danger happening now.

Emergent Mode

If there is emergency happens, your phone will pop up notification.

- To clear signal after you ensure no one is in danger at home, press "Clear".

NOTE: Please contact the user at home then clear signal if appropriate.

- To call the one at home, press "Ask".
- To call EMS immediately, press "911".

